

Quad cane with large or small base



The rotating base swivels and locks into position to accommodate right-handed or left-handed users and is available in small or large for greater stability. Push-button height settings ensure a proper fit, which reduces stress on your wrist, shoulder and back.

- Height adjustable from 28" to 39" (71 to 99 cm)
- 11 push-button height adjustments fit most people 5'-6'5" (152 – 198 cm)*
- Select the push-button setting that puts the cane handle at wrist level
- Shock Absorbing Tips – Ultra-Grip Edge™ with flared rubber tips for extra stability and traction
- K-shaped Quad Base – Small or large base provides greater stability, easily rotates for right or left hand use
- Free Standing – Conveniently stands on its own
- Maximum Capacity – 300 lbs (136 kg)

Walking With a Quad Cane

To walk with a quad cane, follow these simple steps.

1. Hold the cane in one hand. If one leg is weak, hold the cane in the hand opposite the weak leg.
2. Advance the quad cane forward about one arm's length. Be sure all four legs of the quad cane contact the floor to prevent tipping.
3. Step forward with the weak leg.
4. Gently press down into the handle of the quad cane with your hand to help with stability. Advance your other leg just slightly past to first foot.
5. Repeat this cycle.

Be sure your PT helps you set your quad cane up and avoid common mistakes that some folks make when walking with a cane.