

Overbed Table



(typical 3-in-1 commode)

An overbed table allows a stable surface for a variety of uses for users that need to spend time in bed. The table's U-shaped design, height adjustability, and casters allows the table to be rolled under the bed frame and therefore positioned over the user instead of at the side of the bed. This enables easier eating, reading, playing games, writing, electronics use, and more.

Overbed Table Instructions

1. To raise the table to desired height, simply lift up under the table (some tables may require squeezing a tab or lever on the table).
2. To lower the table, squeeze the tab or lever on the table and push down on the table column to the desired height. Do not push down on the tabletop end without the column.
3. The wheels can be locked into place once the table is positioned by depressing the small tab locks on the wheels.